

ERGONEERS

FROM SCIENCE TO INNOVATION

QUICK START D-LAB

Eye Tracking (05.2018)

D-LAB UPDATES

1. Backup your data: http://www.ergoneers.com/faq/index.php?solution_id=1185
2. Download the latest D-Lab setup and install it on your computer, in case you have a wireless system also download D-Lab Mini and install it on the tablet
<http://www.ergoneers.com/faq/index.php?action=news&newsid=3&newslang=en>
3. Please note the Release Notes, which include known & fixed bugs and list new features
<http://www.ergoneers.com/faq/index.php?action=news&newsid=2&newslang=en>


TROUBLE SHOOTING

You can find answers to your questions, as well as instructions, video tutorials and troubleshooting information on the Ergoneers FAQ (<http://www.ergoneers.com/faq>)

In case you experience any further issues, please send us an e-mail to support@ergoneers.com

WHAT IS D-LAB?

- Acquire, synchronize & combine data...
- from multiple sources...
- to generate insights into human behavior....
- While being in a specific situation/solving a task/interacting with a device or using social media.






D-LAB

SIMPLIFIES THE EXPERIMENTAL WORKFLOW

PLAN

- Organise subjects
- Create experimental design
- Configure, calibrate and manage sensors
- Control Marker position




MEASURE

- Observe data acquisition
- Observe and interact with subject
- Comment events
- Trigger experimental conditions to structure long recording




ANALYZE

- Check data quality
- Post-process data
- Visualising data
- Calculate quantitative gaze data to defined AOIs
- Extract data out of immense data collection and create new channels
- Export data to further analyse it statistically





EYE TRACKING WITH DIKABLIS



DIKABLIS GLASSES 3 (2018)



Headunit



eBox

DIKABLIS PROFESSIONAL (2014)

1 Scene camera

- 1920 x 1080 Px (Full HD)
- Opening angle up to 90°
- Adjustable



Nose support

Adjustable for a comfortable fit




2 Eye cameras

- Adjustable angle
- Eye illuminating IR Cameras



EYE TRACKING - WORKING PRINCIPLE

Experimental situation



3 Videos form the Dikablis Eye Tracker in D-Lab



Scene Cam Video




R-Eye video



L-Eye video

EYE TRACKING - WORKING PRINCIPLE




4 point calibration



D-Lab

Pupil detection algorithm tracks the center of the pupil




Video of visual Scene with gaze point



TERMS & DEFINITIONS ISO 15007-1

HOW DO WE DESCRIBE VISUAL BEHAVIOUR?



TERMS & DEFINITIONS ISO 15007-1

HOW DO WE DESCRIBE VISUAL BEHAVIOUR?

Fixation

Alignment of the eyes so that the image of the fixated area of interest falls on the fovea for a given time period (duration from 100ms - 2000ms)



Saccade

Brief, fast movement of the eyes that changes the point of fixation

Areas of Interest (AOI)

Pre-determined areas within the visual scene...any shape



Glance duration

Maintaining of visual gaze within an AOI → Sum over one or more fixations and saccades to one AOI.



TERMS & DEFINITIONS ISO 15007-1

WHAT IS A GLANCE?






D-LAB EYE TRACKING WORKFLOW

PLAN

- Control Marker position
- Adjust Dikablis cameras
- Calibrate Eye Tracker



Recording Devices		
Name		Status
Audio		
Mikrofonarray (Realtek High Definition Audio)	<input checked="" type="checkbox"/>	Connected
Eye Tracker	<input checked="" type="checkbox"/>	
Dikablis Professional	<input checked="" type="checkbox"/>	Connected
Eye Cam - Left		
Eye Cam - Right		

MEASURE

- Monitor quality of calibration
- Correct calibration (optional)



ANALYZE


- Check pupil detection rate
- Post-process data (Marker detector)
- Adjust calibration (optional)
- Draw AOIs
- Calculate quantitative gaze data to defined AOIs
- Calculate statistics (glance duration, number of glances....)



EXCURSION: TASKS & MARKER

EXPERIMENTAL PLAN

ISO/TS 15007-2



Task Definitions	
Name	Network Name
TomTom	TomTom
Navigation	TomTom\Navigation
Town	TomTom\Navigation\Town
Street	TomTom\Navigation\Street
House Number	TomTom\Navigation\House Number

MARKER TECHNOLOGY



TO DRAW DYNAMIC AOIS & AUTOMATIC GLANCE CALCULATION

What does Markers do?


This 16 QR-Codes can be detected by D-Lab and are used as reference points in the experimental environment, to deal with head movements

Dos

- Use each Marker only once in experimental environment
- Attach Markers in same plane (distance from Scene Cam) like the object is, you want to analyse glances to
- Attach Marker fixed on an object, if object will be moved in visual scene
- Print Markers in different sizes (optional)
- Use Offline Marker Detector in Analyse Mode in case of bad online Marker detection
- Marker need enough light to be visible in Scene Camera




MARKER TECHNOLOGY: EXAMPLE



Head movement to the left




Marker Oslo moves to right in Scene Cam image like the "W" button (framed by the purple AOI)



DIKABLIS CALIBRATION – WORKFLOW


D-LAB PLAN MODE



EYE TRACKING DATA ANALYSIS – WORKFLOW

D-LAB ANALYSE MODE


1. Check Pupil Detection Rate



For rates < 80 % use “Automatic Pupil Detection” offline  in Eye Tracker visualization or adjust center of pupil manually (labor & time intensive).



2. Recalibration (Opt.)



ONLY for advanced users:
It is the option to adjust Distance (+/-), Position (→/←) and Rotation of Eye Video in relation to Scene Video



3. Offline Marker Detection (Opt.)



Data Analysis Tab → Marker Detection
Fast Detection
Normal Detection
Exhaustive detection (time intensive)



6. Calculate Eye Tracking Statistics

Data Analysis Tab → Eye Tracking Statistics



Start with from left to right and select Subject/Groups, Time interval (Tasks), AOIs and Metrics that like to be calculated for the selected combinations.




7. Export Data



5. Calculate Glances

Data Analysis Tab → AOI Management


Use all four buttons from left to right to get the most valid eye tracking data according to ISO 15007-1.



4. Define AOIs

Eye Tracker Visualization

Navigate to a frame that shows the object of interest & a recognized marker to bound it to



DIKABLIS METRICS & DATA CHANNELS

EYE TRACKING METRICS

$$\text{AOI attention ratio [\%]} = \frac{\text{Total Glance Time on AOI}}{\text{Duration of Task}} * 100$$

Percentage of time of the selected time interval, in that glances are within an AOI (or set of related AOIs)

Maximum/Minimum Glance Duration [s]:

Longest/Shortest glance duration at a specific AOI in the selected time interval (task/subtask/sub-subtask)

$$\text{Glance location probability [\%]} = \frac{\text{Number of glance to at AOI}}{\sum \text{Number of glances to AOI1, AOI2, ... , AOIn}} * 100$$

Probability that the test person looks at an AOI (or AOI set) during a particular time interval.

$$\text{Link value probability [\%]} = \frac{\text{Number of transitions from A to B} + \text{Number of transitions from B to A}}{\text{Number of all transitions in given time}} * 100$$

Probability of glance transitions between two different AOIs (A&B). To calculate this metric, two AOIs have to be selected for calculation.

$$\text{Percentage transition time [\%]} = \frac{\sum \text{Duration 1 outside all selected AOIs, duration 2 outside all selected AOIs, ..., duration n}}{\text{Duration of selected time interval}} * 100$$

Represents the percentage of time (100% is the duration of selected time interval), when the gaze is not in any selected AOI.

EYE TRACKING METRICS

$$\text{Horizontal/Vertical Eye Activity [Pixel]} = \sqrt{\frac{\sum(x-\bar{x})^2}{n-1}}$$

Standard deviation of the X or Y pupil coordinate in the coordinate system of the scene image. This metric is AOI independent and indicates the visual activity during the selected time interval (task).

$$\text{PERCLOS left/average/right [\%]} = \frac{\text{number of frames without detected pupil}}{\text{total number of frames in time interval}} * 100$$

The Percentage of Eye Closure is calculated as number of frames with no detected pupil, divided by the total number of frames in the selected time interval.

Mean fixation duration left/right [ms]

Length of time that a glance is fixed on a particular AOI in seconds (such fixations are calculated according to the principle of Dario D. Salvucci and Joseph H. Goldberg, see also: <http://doi.acm.org/10.1145/355017.355028>)

Number of fixations

Number of fixations on the AOI for the selected time interval.

Mean saccade duration left/right [ms]

Sum of saccade durations in selected time interval divided by number of saccades in selected time interval.

EYE TRACKING METRICS

Mean saccade angle left/right [deg]

Sum of saccade angles in selected time interval divided by number of saccades in selected time interval (assumption: 1deg = 2,5px movement of crosshair in Scene Image)

Number of saccades right/left

Number of saccades for the selected time interval.

Time to first glance [s]

Duration from the beginning of the selected time interval until the first glance to an AOI.

Glance duration Percentile [s]

Enter a percentile you are interested into the empty field and D-Lab calculated the percentiles glance duration. E.g. type 0.95 and D-Lab calculates the duration in seconds, that 95% of all other measures glance durations fall short of this duration.

EYE-TRACKING RAW DATA – EYE-DATA

The information provided under “Eye-Data” is the same for every eye (left / right). It is shown in the picture below.

Left Eye	
Pupil X	Real
Pupil Y	Real
Pupil Area	Real
Pupil Width	Real
Pupil Height	Real
Index Of Cognitive Activity Average	Real
Index Of Cognitive Activity Sample	Real
Saccades	
Saccades	Real
Saccades Duration	Real
Saccades Angle	Real
Fixations	
Fixations	Real
Fixations Duration	Real

Saccades and fixations are calculated for every eye separately. The calculation is based on the criteria set in the options menu (File → Options):

Saccade / Velocity Threshold

→ sets the limit for the movement speed of the pupil to differentiate between a saccade and a fixation.

Movement speed > Threshold: interpretation as saccade

Movement speed < Threshold: interpretation as fixation.

Pixels per degree Coefficient

→ As the distance of the eye camera to the eye can vary, the image size of the eye can be different. Therefore it is necessary to adjust the detected eye-movement in Pixel to an appropriate movement in degree.

EYE-TRACKING RAW DATA – EYE-DATA

The information provided under “Eye-Data” is the same for every eye (left / right). It is shown in the picture below.

Left Eye

Pupil X	Real
Pupil Y	Real
Pupil Area	Real
Pupil Width	Real
Pupil Height	Real
Index Of Cognitive Activity Average	Real
Index Of Cognitive Activity Sample	Real

Saccades

Saccades	Real
Saccades Duration	Real
Saccades Angle	Real

Fixations

Fixations	Real
Fixations Duration	Real

0 or 1, depending if there is a saccade detected in the current frame. (Boolean)

Duration of the Saccade in [s].

Angle of the saccade in [°].

0 or 1, depending, if there is a fixation detected in the current frame. (Boolean)

Duration of the Fixation in [s].

EYE-TRACKING RAW DATA – FIELD DATA

- Field Data

- Scene Cam

- Gaze

Gaze X	Real
Gaze Y	Real

X / Y Position of the gaze (crosshair) in the coordinate system of the Scene Camera is given in [px]. Gaze calculation is also possible for gaze positions beyond the scene image e.g. a gaze position **G1(-25/1010)**

- Markers


- Brasilia

Brasilia X1	Real
Brasilia Y1	Real
Brasilia X2	Real
Brasilia Y2	Real
Brasilia X3	Real
Brasilia Y3	Real
Brasilia X4	Real
Brasilia Y4	Real

- Gaze in Marker Coordinates

- Brasilia

Pupil X	Real
Pupil Y	Real





EYE-TRACKING RAW DATA – FIELD DATA

- Field Data
 - Scene Cam
 - Gaze
 - Gaze X
 - Gaze Y
 - Markers
 - Brasilia
 - Brasilia X1
 - Brasilia Y1
 - Brasilia X2
 - Brasilia Y2
 - Brasilia X3
 - Brasilia Y3
 - Brasilia X4
 - Brasilia Y4
 - Gaze in Marker Coordinates
 - Brasilia
 - Pupil X
 - Pupil Y

Detected Marker: Markers need to be detected by D-Lab at least once for every recording, before they show up in the treeview.


X / Y Position of the 4 corners of the detected marker in the coordinate system of the Scene-Camera, given in [px]. The Markerorigin (P1) is marked with a grey dot.




EYE-TRACKING RAW DATA – FIELD DATA

The unit of the coordinate system and thus the unit in which the position of the Gaze in Marker-coordinates is given, is the edge length of the marker as it is detected in the scene image. E.g. a Gaze Position of $X = 1.5$ and $Y = 2$ means, that the crosshair is located 1.5x the size of the marker in x-direction, and 2x the size of the marker in y-direction. Thus, even when the Scene Camera moves closer to the marker, the relation within the coordinate system stays the same.

- Field Data
 - Scene Cam
 - Gaze
 - Gaze X
 - Gaze Y
 - Markers
 - Brasilia
 - Brasilia X1
 - Brasilia Y1
 - Brasilia X2
 - Brasilia Y2
 - Brasilia X3
 - Brasilia Y3
 - Brasilia X4
 - Brasilia Y4
- Gaze in Marker Coordinates
 - Brasilia
 - Pupil X
 - Pupil Y



X / Y Position of the Gaze-Crosshair in the coordinate system of the respective marker (e.g. Brasilia). The unit is the detected edge-length of the marker [-].



Gaze in Marker Coordinates (-1/1,5)
 Gaze in Scene Cam Image (500/300)

MORE MANUALS AVAILABLE HERE:

<http://www.ergoneers.com/faq/index.php?action=artikel&cat=7&id=52&artlang=en>