

Produkt: Hörfitness® App



Verfügbare Plattformen: Android | iOS
Verfügbare Hardware: Tablet (empfohlen!)
Smartphone

Mindestanforderungen Tablet:

Android-Version: 10 oder höher
Display-Größe: 10 Zoll oder größer (empfohlen)
Display-Auflösung: 1920 x 1080 dpi (Full-HD)
Hauptspeicher (RAM): 3 GB oder größer
Speicher (HD): 64 GB oder größer

Aktuelle Informationen zu Versionsständen und technischen Anforderungen unter
<https://www.meditech.de/software-downloads>

MediTECH-Infokanal (öffentlich):
<https://www.meditech.de/slides/horfitness-infokanal-129>

MediTECH-Kundenkanal (exklusiv für Kaufkunden):
<https://www.meditech.de/slides/horfitness-kundenkanal-121>

Product: **Audiofitness app**

Available platforms: Android | iOS
Available hardware: Tablet (recommended!)
Smartphone



Minimum tablet requirements:

Android version: 10 or higher
Display size: 10 inch or larger (recommended)
Display resolution: 1920 x 1080 dpi (Full-HD)
Main memory (RAM): 3 GB or larger
Memory (HD): 64 GB or larger

Current information on version statuses and technical requirements can be found at https://www.meditech.de/en_EN/software-downloads

MediTECH info channel (public):

https://www.meditech.de/en_US/slides/audiofitness-info-channel-129

MediTECH customer channel (exclusive for Purchasing customers):

https://www.meditech.de/en_US/slides/audiofitness-customer-channel-121