

**Produkt:** **Body & Mind App**

Verfügbare Plattformen: Android | iOS

Verfügbare Hardware: Tablet (empfohlen!)  
Smartphone



**Mindestanforderungen Tablet:**

Android-Version: 10 oder höher

Display-Größe: 10 Zoll oder größer (empfohlen)

Display-Auflösung: 1920 x 1080 dpi (Full-HD)

Hauptspeicher (RAM): 3 GB oder größer

Speicher (HD): 64 GB oder größer

Aktuelle Informationen zu Versionsständen und technischen Anforderungen unter

<https://www.meditech.de/software-downloads>

MediTECH-Infokanal (öffentlich):

<https://www.meditech.de/slides/body-mind-app-infokanal-67>

**Product:** **Body & Mind App**

Available platforms: Android | iOS

Available hardware: Tablet (recommended!)  
Smartphone



**Minimum tablet requirements:**

Android version: 10 or higher  
Display size: 10 inch or larger (recommended)  
Display resolution: 1920 x 1080 dpi (Full-HD)  
Main memory (RAM): 3 GB or larger  
Memory (HD): 64 GB or larger

Current information on version statuses and technical requirements can be found at [https://www.meditech.de/en\\_EN/software-downloads](https://www.meditech.de/en_EN/software-downloads)

MediTECH info channel (public):

[https://www.meditech.de/en\\_EN/slides/body-mind-app-infokanal-67](https://www.meditech.de/en_EN/slides/body-mind-app-infokanal-67)